



## Canespro Fungal Nail Treatment Set

### Uses

Removes the infected parts of the nail in 2-3 weeks.

### Dosage

Step 1: Apply a thin layer of the ointment to the infected nail.

Step 2: Cover the treated nail with the plaster provided and keep covered for 24 hours.

Step 3: Remove the plaster and soak the nail in warm water for approximately 10 minutes. Gently remove the softened parts of the nail using the plastic scraper.

Repeat steps 1, 2 and 3 daily for 2-3 weeks.

### Format and Ingredients

Treatment set containing:

10g tube Canespro Urea 40% ointment, 22 waterproof

A division of  
 GROUPE  
EUROCOM

 PAGB

 **Communications  
International  
Group**